SAM® IO Driver | Instructions for Use

PRODUCT INFORMATION

DESCRIPTION:

The SAM IO Intraosseous Access System utilizes a manually operated driver handle. The driver functions by continuously actuating (repeatedly compressing) the trigger assembly. Repeated, full trigger actuation creates rotational spin with a passive ratcheting mechanism designed to place the SAM IO needle assembly.

PRODUCT INFORMATION:

- Driver part number: IO700-EN
- Applied Parts: SAM IO Intraosseous Access needles; 15 mm, 25 mm, 45mm

SAFETY INFORMATION:

- Indications, contraindications, warnings, precautions, and other safety information are contained in the Instructions for Use for the SAM IO Intraosseous Access System.
- Please consult the Instructions for Use for the SAM IO Intraosseous Access System before applying. If there are questions, immediately contact your local SAM Medical sales representative.
- Additional product information can be found at sammedical.com

IMPORTANT INFORMATION FOR USERS:

For SAM IO Intraosseous Access System products to perform properly, the following conditions are recommended.

- Use this product only in accordance with this manual and applicable product labeling.
- Adjustments, modifications, technical maintenance or repairs are not allowed.
- Do not connect this product or its components to products not recommended by SAM Medical.
- Visually inspect driver for deficiencies before use. If driver is soiled, discard or clean according to the method below:
 - Maintain personal protective equipment precautions.
 - Wipe exterior surfaces of the SAM IO driver with soft, clean moistened cloth.
 - Use a soft bristled brush dampened with a 1% Alconox solution to remove visible soil. Clean all surfaces that make contact with the hand during use, with particular attention paid to the rubberized grip on the front of the driver.
 - Rinse the IO driver under running tap water, but do not rinse the driver's internal components.
 - Dry driver with a soft, clean cloth.

IMPORTANT INFORMATION FOR USERS CONTINUED:

- After cleaning, inspect to ensure no visible debris remains and no damage has occurred to the driver. If soils remain on or inside the driver, repeat steps above or follow the additional method below:
 - Submerge the driver in 1% Alconox solution for five minutes.
 - Use a soft bristled brush dampened with a 1% Alconox solution to remove any remaining visible soil.
 - Rinse the driver's exterior surfaces and internal components with tap water, ensuring the driver is flushed of 1% Alconox solution.
 - Dry driver with a soft, clean cloth.
 - If soils remain on or inside driver, discard the SAM IO Driver.
- Spray or wipe exterior surface with anti-microbial solution or towelettes (e.g. Super Sani-Cloth® or equivalent intermediate level EPA registered disinfectant agent) following the manufacturer's recommendations.
- Actuate the driver 1-3 times and ensure no tactile or visible damage is felt or seen before returning the driver to its storage location. If tactile or visible damage is noted, discard the driver.
- Unless mechanically damaged, there is no known expiration of the driver. Mechanical damage ending the useful life of the product will be readily identified as fractures, bends, or other visibly obvious damage. Before attachment of the needle, actuate driver 1-3 times and if no tactile or visible damage is felt or seen the driver is viable for use.


Do not use excessive force during insertion. The mechanical rotation of the needle by handle actuation and the cutting edge of the needle should be the PRIMARY mechanisms to penetrate the bone, NOT the downward pressure.

If driver is unavailable, needle assembly can be used manually. Please refer to SAM IO Needle Instructions for Use (step 5).

Not for Sternal use.

SYMBOL GLOSSARY:

For prescription use only.	
Indicates the need for the user to consult the instructions for use.	
Indicates the medical device manufacturer, as defined in EU Directive 93/42/EEC.	
Indicates the manufacturer's batch code so that the batch or lot can be identified.	
Indicates the manufacturer's catalog number so that the medical device can be identified.	
Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	
Not for Sternal use.	
Indicates a medical device that has not been subjected to a sterilization process.	

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