



SAM
MEDICAL

SAM[®] IO

Intraosseous Access System

**MULTI-USE.
DURABLE.
LIGHTWEIGHT.**



PROVEN TO SURVIVE

Engineered for rapid application, SAM IO's innovative & lightweight design is battery free, and enables easier, faster training and intuitive use. Meets MIL-STD 810G for durability.

ENGINEERED FOR SURVIVAL

MANUALLY OPERATED DRIVER

SAM IO is a manually operated intraosseous access system. Catheter placement is achieved by continuously actuating (repeatedly compressing) the driver's trigger assembly while gently guiding needle assembly into position. Repeated, full trigger actuation creates rotational spin of the needle assembly which, when combined with gentle downward pressure, results in controlled IO placement. Once needle assembly is properly positioned, the stylet is removed to expose standard Luer-lock for extension set connection. With the extension set connected, aspiration verification, flushing and selected treatment(s) may commence.





SAM[®] IO

SYSTEM OVERVIEW

- 1 SAM IO NEEDLES**
A stylet and catheter assembly secured by a press-fit connection. 15 mm, 25 mm, and 45 mm lengths are intended to address a broad range of patient demographics and multiple application sites.
- 2 NEEDLEWISE[®] BLOCK**
Included in every Needle Set. Single-use sharps containment device for the stylet.
- 3 SAM IO DRIVER**
Light weight (2.5 oz), a multi-use manually operated intraosseous access driver. The driver, made of a rugged Polycarbonate, was tested for over 10,000 actuations and meets MIL STD 810G.

- 4 SAM STABILIZER**
Cylinder design intended to protect the catheter from dislodgement.
- 5 EXTENSION TUBING**
Connects to the catheter for flexible infusion.
- 6 SAM IO NEEDLE ADAPTOR**
Attaches to the SAM IO Driver to accommodate alternative IO needles.



PRODUCT SPECIFICATIONS

DRIVER SIZE & WEIGHT*

Length: 5.5"
Width: 1"
Height: 4.25"
Weight: 2.5 oz

MODELS AVAILABLE

Driver:
IO700-EN

Stabilizer:
1-Pack IO740-1P-EN
5-Pack IO740-5P-EN

Adaptor:
IO742-EN

Needles:

15 mm 1-Pack IO705-1P-EN
15 mm 5-Pack IO705-5P-EN
25 mm 1-Pack IO706-1P-EN
25 mm 5-Pack IO706-5P-EN
45 mm 1-Pack IO707-1P-EN
45 mm 5-Pack IO707-5P-EN

GLOBAL SALES

U.S. Toll Free: 800.818.4726
U.S. Local: 503.639.5474
U.S. Fax: 503.639.5425
sammedical.com
sammedical.store




SAM Medical Products[®]
12200 SW Tualatin Rd., Ste 200
Tualatin, OR 97062 USA

* Size and weight vary by configuration and manufacturing process.
SAM[®] is a registered trademark in the United States and other countries.
©2020 SAM Medical Products. All rights reserved.